

# Detergents/Floc Testing

## DETERGENTS

### LaMotte Drop Count Kit

Easy to use drop count kit

**Measuring range:** Each drop equivalent to 1.0ppm

Supplied with reagents for approx. 60 tests

**Reference:** LE/6605-0217                      £115.00



LaMotte Drop Count Kit

## FLOC TESTING

### Phipps and Bird Standard Jar/Floc Tester

The Phipps and Bird unit is available in a four or six place format. Constructed on a powder coated steel uniframe chassis the stainless steel paddles are spaced 150mm apart and are adjustable to a maximum depth of 230mm.

#### Features:

- Stirring speed 1 to 300rpm.
- Illuminated base.
- Digital readout

Each unit is supplied with a dust cover and either 1 litre standard glass beakers or 2 litre square acrylic testing jars. The four position unit also has the option of a heavy duty carrying case.



Phipps and Bird Standard Jar/Floc Tester

- Six place unit with 6 round beakers.  
**Reference:** PB/7790-702                      £3010.00
- Six place unit with 6 square beakers.  
**Reference:** PB/7790-712                      £3250.00
- Four place unit with 4 round beakers.  
**Reference:** PB/7790-752                      £2760.00
- Four place unit with 4 square beakers.  
**Reference:** PB/7790-762                      £3010.00
- Carrying case for four place unit.  
**Reference:** PB/7791-020                      £550.00



FOR DISSOLVED OXYGEN ANALYSIS GO TO PAGE 96