

# Fluoride/Formaldehyde

## FLUORIDE

### Camlab CW2080 Colorimeter

Single parameter handheld colorimeter. Microprocessor controlled for accurate and reproducible results.

**Measuring range:** 0.2 – 0.5mg/l F

**Wavelength:** 528nm

**Features:**

- Automatic zero calibration
- LED light source with long term stability
- Large digital display
- Auto switch off

**Reference:** CW/20.93.00

£398.00



Camlab CW2080 Colorimeter

### Hach Pocket Colorimeter

Single parameter handheld colorimeter. Fixed wavelength and three key easy to operate design.

**Measuring range:** 0. – 2.0mg/l F<sup>-</sup> in steps of 0.1

**Chemistry:** SPADNS

**Features:**

- Silicon photodiode
- LED light source
- IP67 casing
- 2000 test battery life

Supplied with reagents for 100 tests, accessories, carrying case and instructions

**Reference:** HH/58700-05

£362.30



Hach Pocket Colorimeter

## FORMALDEHYDE

### Hach Colour Match Kit Model FM-2

Easy to use colour match kit complete with reagents for 100 tests, instructions and carrying case.

**Increments:** 0, 0.5, 1 and 5mg/l as CH<sub>2</sub>O

**Chemistry:** MBTH

**Reference:** HH/22672-00

£64.50



Hach Colour Match Kit Model FM-2