

Hazen (Colour)/Hydrazine

HAZEN (COLOUR)

Camlab CW2100 Colorimeter

Single parameter handheld colorimeter. Microprocessor controlled for accurate and reproducible results.

Measuring range: 0 – 500 Pt-Co units

Wavelength: 470nm

Features:

- Automatic zero calibration
- LED light source with long term stability
- Large digital display and auto switch off

Reference: CW/20.72.40 £270.00

Colour Standard Solutions

All supplied in 1 litre bottles.

0 Hazen Units	
Reference: RN/HAZ0	£49.17
10 Hazen Units	
Reference: RN/HAZ10	£49.17
25 Hazen Units	
Reference: RN/HAZ25	£49.17
40 Hazen Units	
Reference: RN/HAZ40	£56.35
50 Hazen Units	
Reference: RN/HAZ50	£56.35
80 Hazen Units	
Reference: RN/HAZ80	£61.49
100 Hazen Units	
Reference: RN/HAZ100	£61.49
250 Hazen Units	
Reference: RN/HAZ250	£61.49
500 Hazen Units	
Reference: RN/HAZ500	£73.27



Camlab CW2000 Colorimeter

HYDRAZINE

Camlab CW2110 Colorimeter

Single parameter handheld colorimeter. Microprocessor controlled for accurate and reproducible results.

Measuring range: 0.05 – 0.5mg/l as N₂H₄

Wavelength: 470nm

Features:

- Automatic zero calibration
- LED light source with long term stability
- Large digital display and auto switch off

Reference: CW/20.73.00 £270.00



Camlab CW2000 Colorimeter