

Hach Colour Comparator Kit Model IR-24

Colour disc kit with reagents for 100 tests and instructions. For low range soluble iron.

Measuring ranges:

0 – 0.2mg/l in steps of 0.002 dissolved Fe²⁺/³⁺
0 – 1.0mg/l in steps of 0.01 dissolved Fe²⁺/³⁺

Chemistry: Ferrozine

Reference: HH/02556-00 £98.50



Hach Colour Comparator Kit

Hach Colour Comparator Kit Model IR-21

Colour disc kit with reagents for 100 tests and instructions. For low range soluble iron using the TPTZ method.

Measuring ranges:

0 – 0.1mg/l in steps of 0.01 dissolved Fe²⁺/³⁺
0 – 1.2mg/l in steps of 0.05 dissolved Fe²⁺/³⁺

Chemistry: TPTZ

Reference: HH/22993-00 £95.10

Hach Cube Kit

Simple to use colour comparison kit complete with reagents for 50 tests and instructions.

Measuring range: 0 – 5mg/l in steps of 1 dissolved Fe²⁺/³⁺

Chemistry: 1,10 Phenanthroline

Reference: HH/14008-00 £22.90



Hach Cube Kit

Hach Cube Kit

Simple to use colour comparison kit complete with reagents for 50 tests and instructions.

Measuring range: 0 – 10mg/l in steps of 2 dissolved Fe²⁺/³⁺

Chemistry: 1,10 Phenanthroline


Reference: HH/25435-00 £28.50

Quantofix Test Papers for Iron

Easy to use test papers supplied in boxes of 100.

Graduations: 0-5-20-50-100-250-500-1000mg/l Fe²⁺/³⁺

Reference: MN/913.02 £14.30



Quantofix Test Papers

Quantofix Test Papers for Iron

Easy to use test papers supplied in boxes of 100.

Graduations: 0-2-5-10-25-100mg/l Fe²⁺/³⁺

Reference: MN/913.08 £14.30