

MANGANESE

Camlab CW2140 Colorimeter

Single parameter handheld colorimeter. Microprocessor controlled for accurate and reproducible results of manganese.

Measuring range: 0.05 – 4mg/l Mn²⁺ /³⁺

Wavelength: 430nm

Features:

- Automatic zero calibration
- LED light source with long term stability
- Large digital display
- Auto switch off

Reference: CW/20.94.00

£270.00



Camlab CW2000 Colorimeter

Hach Pocket Colorimeter

Single parameter handheld colorimeter. Fixed wavelength and three key easy to operate design.

Measuring range: 0. – 20.00mg/l in steps of 0.1 as Mn

Chemistry: Cold Periodate

Features:

- Silicon photodiode
- LED light source
- IP67 casing
- 2000 test battery life

Supplied with reagents for 100 tests, carrying case, accessories and instructions

Reference: HH/58700-15

£362.30



Hach Pocket Colorimeter

Hach Pocket Colorimeter

Single parameter handheld colorimeter. Fixed wavelength and three key easy to operate design.

Measuring range: 0 – 0.7mg/l in steps of 0.05 as Mn

Chemistry: PAN

Features:

- Silicon photodiode
- LED light source
- IP67 casing
- 2000 test battery life

Supplied with reagents for 100 tests, carrying case, accessories and instructions

Reference: HH/58700-18

£362.30



Hach Pocket Colorimeter