

Lovibond Comparator Kit

Checkit comparator kit complete with colour disc, reagents and carrying case.

Measuring range: 0 – 0.5mg/l N

Reference: TT/14.73.00 £29.50

Hach Colour Comparator Kit Model NI-6

Colour disc kit with reagents for 100 tests and instructions.

Measuring ranges: 0 – 100mg/l in steps of 2 as NO_2^-
0 – 20000mg/l in steps of 40 as NO_2^-

Chemistry: Ferrous Sulphate

Reference: HH/02240-00 £110.30

Hach Colour Comparator Kit Model NI-15

Colour disc kit with reagents for 100 tests and instructions.

Measuring ranges: 0 – 0.5mg/l in steps of 0.01 as NO_2^- -N

Chemistry: Diazotization

Reference: HH/21820-00 £72.10

Hach Cube Kit

Simple to use colour comparison kit complete with reagents for 50 tests and instructions.

Measuring range: 0 – 1.0mg/l in steps of 0.2 as NO_2^- -N

Chemistry: Diazotization

Reference: HH/20596-00 £24.30

Quantofix Test Papers for Nitrite

Easy to use test papers supplied in boxes of 100.

Graduations: 0-1-5-10-20-40-80 mg/l as NO_2^-

Reference: MN/913.11 £16.80

Quantofix Test Papers for Nitrite

Easy to use test papers supplied in boxes of 100.

Graduations: 0-0.1-0.3-0.6-1-2-3g/l as NO_2^-

Reference: MN/913.22 £16.80



Lovibond Comparator Kit



Hach Colour Comparator Kit



Quantofix Test Papers