

Autoclaves



Swiftlock Compact Autoclaves
23, 35 and 50 litre

AUTOCLAVES

Swiftlock Compact Autoclaves 23, 35 and 50 litre

The new ranges of autoclaves offers guaranteed precision and versatility for a wide range of applications and are designed to meet the demanding requirements of current sterilisation procedures.

Features:

- 4 program micro-control
- 316L stainless steel pressure vessels
- Temperature range 100 to 138°C (0.2-2.4 bar)
- Ultra-quick new door lock mechanism
- Conforms to national and international standards
- Automatic water fill
- Built in drying system
- Water recycling as standard

The Swiftlock automatic water fill/recirculation system together with extremely versatile programming capability provides exceptional sophistication in a small autoclave. The integral reservoir, once filled, provides sufficient water for up to 20 cycles; moreover, the same container also acts as collection vessel for both condensate and water (which is returned from the autoclave chamber following the sterilisation process). Thus the system is completely self-contained and connection to a drain or waste container is unnecessary.

SPECIFICATIONS

Capacity (litres)	Loading	Chamber size (dia x depth mm)	Reference	Price
23	Front	250 x 430	AL/AMB220	£4180.00
40	Front	350 x 400	AL/AMB230	£5320.00
60	Front	350 x 600	AL/AMB240	£6635.00
60	Top	350 x 570	AL/AMA240	£6440.00

Optional accessories

Integral data printer

Reference: AL/AAR130 £780.00

Load sensor process timing

Reference: AL/AAR014 £245.00



Acculab Pocket-Pro Balance

BALANCES

Acculab Pocket-Pro Balance

A battery operated portable balance, ideal for sample preparation, test kits and special applications. Has a rugged ABS housing with a slide closure mechanism that locks the weighing pan during transport as well as protecting the backlit display. Also includes low battery, overload and underload indicators, weighing units in grams and ounces, simple 3 key operation. Weighing tray, battery and calibration weight included. Capacity 200g x 0.1g

Reference: SR/PP-201

£40.00