


Balances



Sartorius Moisture Balance

Sartorius Moisture Balance

The MA35 is an easy to use moisture balance with automatic end point determination. Is idea for routine analysis where accuracy and ease of use are important.

- Easy to understand icon lead menu.
- Two metal tubular heating elements providing 360 watts of power.
- Automatic end point determination.
- Weighing capacity 35g, accuracy 1mg.
- Temperature range 40 to 160°C.
- RS-232C interface.
- Built in weighing system to provide measurement accuracy.
- Display of results in % moisture, % dry weight, % RATIO and g residue.

Reference	Price
SR/MA35	£1295.00



Denver Moisture Balance

Denver Moisture Balance

The IR-35 is suitable for all applications and offers great performance especially matched to the daily needs of the laboratory. Its easy operational concept provides a large and clearly arranged display that gives all important information at a glance.

- Automatic end point determination.
- Weighing capacity 35g
- Accuracy 1mg
- Temperature range 40 to 160°C
- RS-232C interface
- Fully automatic analysis mode
- Three easy steps from taring to analysis
- Display of analysis in moisture, dry weight, RATIO and residual weight

Reference	Price
DR/IR-35	£1225.00