

## Universal Water Baths

Designed for direct immersion of samples. Has integrated electronic control. The working area is completely unobstructed and the baths are designed to be placed 'end-on' for convenient access and observation, and for efficient utilisation of bench space.

### Temperature Range

- Ambient +5 to 99°C (stability  $\pm 0.2^\circ\text{C}$ )

### Digital Setting and Display

- Ensure reproducibility of set temperature, for repeatability of procedures
- Constant display of actual liquid temperature
- 2-digit easy-to-read LED display

### Base Trays

An accessory tray is required if flat bottomed flasks are to be placed directly on the base of sub universal baths

### Construction

- Consists of a stainless steel tank in an outer case, with heater and temperature sensors mounted underneath the tank
- Recessed control panel protects against spillage

### Safety

- User-resettable, adjustable over-temperature cutout



### Accessories

A full range of accessories are available for these products:

- Gabled Lids
- Flat Lids with Holes
- Raised Shelves
- Base Trays
- Polypropylene Spheres
- Tube Racks

Please call for more information.

## SPECIFICATIONS

	GT/SUB6	GT/SUB14	GT/SUB28	GT/SUB36
Temperature range:	ambient +5 to 99°C			
Temperature setting range:	15 to 99°C			
Stability (DIN 58966):	$\pm 0.2^\circ\text{C}$			
Temperature setting:	digital			
Temperature display:	two digit LED			
Heater power:	0.35kW	0.7kW	1.4kW	2kW
Overall consumption:	0.35kW	0.7kW	1.4kW	
	2kW			
Tank capacity:	6 litres	14 litres	28 litres	36 litres
Working area, lxwx d:	150x300x140mm	325x300x140mm	505x300x190mm	635x300x190mm
Overall dimensions, lxwx d:	205x325x275mm	380x325x275mm	555x325x300mm	720x325x300mm

\* A lid or layer of polypropylene spheres must be used when operating at 60°C or above.

## ORDERING INFORMATION

Universal water bath, 6 litres	GT/SUB6	£425.00
Universal water bath, 14 litres	GT/SUB14	£468.00
Universal water bath, 28 litres	GT/SUB28	£552.00
Universal water bath, 36 litres	GT/SUB36	£673.00
Base tray to fit 6 litre bath	GT/SBT6	£34.00
Base tray to fit 14 litre bath	GT/SBT14	£35.00
Base tray to fit 28 litre bath	GT/SBT28	£36.00
Base tray to fit 36 litre bath	GT/SBT36	£38.00