

# Weighing Boats/Beakers

## Weighing boats

Sufficiently rigid to hold materials on the pan, yet flexible to enable folding for transfer. Excellent for both powders and liquids, they are available in three sizes in black or white for visible contrast with your sample. Made from antistatic polystyrene for easy dispensing.

Reference	Description	Pack Qty	Price
LSE.250	5ml white	1000	£34.00
LSE.251	5ml black	1000	£34.00
LSE.252	30ml white	500	£25.00
LSE.253	30ml black	500	£25.00
LSE.254	100ml white	250	£26.00
LSE.255	100ml black	250	£26.00




Weighing boats

## BEAKERS

### Camlab Economy Beakers

Borosilicate glass beakers with graduations and spout.

Reference	Description	Pack Qty	Price
EU/A2001/04B	50ml	6	£8.77
EU/A2001/05B	100ml	6	£9.22
EU/A2001/07B	250ml	6	£9.90
EU/A2001/08B	400ml	6	£13.50
EU/A2001/09B	600ml	6	£16.42
EU/A2001/11B	1000ml	6	£28.12




Camlab Economy Beakers

### Duran Glass Beakers

Low form beakers borosilicate glass beakers with spout and graduations. DIN12331, ISO3819.

Reference	Description	Pack Qty	Price
SC/2110607	5ml	10	£27.10
SC/2110608	10ml	10	£30.50
SC/2110614	25ml	10	£26.24
SC/2110617	50ml	10	£25.40
SC/2110624	100ml	10	£25.30
SC/2110636	250ml	10	£27.30
SC/2110641	400ml	10	£32.80
SC/2110648	600ml	10	£40.80
SC/6110653	800ml	10	£56.20
SC/2110654	1000ml	10	£64.60
SC/2110663	2000ml	10	£113.00
SC/2110673	5000ml	1	£33.36
SC/2110686	10000ml	1	£92.18



Duran Glass Beakers