


Bottles/Buckets



Glass Reagent Bottles 'Polystop'

Glass Reagent Bottles 'Polystop'

Conform to ISO published recommendations for shape and size. The angle of the shoulders, combined with the dimensions of the ground necks, ensure that contents are poured without 'gulping', thus contributing to safety standards in the laboratory.

Reference	Capacity size	Stopper	Pack qty	Price
Clear glass				
JP/RP705/A	30ml	16/16	10	£19.00
JP/RP705/B	50ml	16/16	10	£19.40
JP/RP705/C	100ml	16/16	10	£23.00
JP/RP705/D	250ml	24/20	10	£27.00
JP/RP705/E	500ml	24/20	2	£7.00
JP/RP705/F	1000ml	29/32	2	£10.20
Amber glass				
JP/RP705/AM/B	50ml	16/16	10	£36.00
JP/RP705/AM/C	100ml	16/16	10	£40.00
JP/RP705/AM/D	250ml	24/20	10	£40.00
JP/RP705/AM/E	500ml	24/20	5	£13.10
JP/RP705/AM/F	1000ml	29/32	2	£17.10



Wash Bottles

Wash Bottles

Polythene, oval shaped with cap and tube.

Reference	Capacity	Price
PWB.3D	250ml	£2.50
PWB.3F	500ml	£3.10

BUCKETS

Stainless Steel Bucket

Manufactured from 304-grade stainless steel, graduated in 1-litre intervals, complete with handle.



Stainless Steel Bucket

Reference	Capacity	Price
LHW.0200	6 litres	£41.00
LHW.0201	8 litres	£56.00
LHW.0202	10 litres	£58.00
LHW.0203	12 litres	£62.00
LHW.0204	15 litres	£65.00



FOR PH BUFFERS SEE PAGES 110 TO 111