

Calibrex digital bottle top dispensers

20ml, 50ml, 100ml adjustable dispensers

Three sizes offer larger volume dispensing with similar precision and accuracy as the Calibrex 520 models.

Ordering information

Range ml	Increments ml	Accuracy (rel error %)	Precision (CV%)	Adjustable volume	Price
2-20	0.1	<±0.7	<0.1	SX/521.020	£270.00
10-50	1	<±0.7	<0.1	SX/521.050	£300.00
10-100	1	<±0.7	<0.1	SX/521.100	£390.00

All 521 units are supplied with adaptors to fit 32, 38, 40 and 45mm screw threads.

Amber glass bottles, square

Reference	Volume	Neck diameter	Price
SX/314-100	100ml	32mm	£5.60
SX/314-250	250ml	32mm	£4.60
SX/314-500	500ml	32mm	£5.80
SX/314-1000	1000ml	45mm	£7.90
SX/314-2500	2500ml	45mm	£10.50

Polyethylene bottles, square

SX/315-250	250ml	25mm	£2.00
SX/315-500	500ml	25mm	£3.20
SX/315-1000	1000ml	25mm	£6.50
SX/315-2500	2500ml	45mm	£7.70

Adaptors, polypropylene, screw neck

SX/GLP22	fits 22mm ø screw necks	£5.50
SX/GLP25	fits 25mm ø screw necks	£5.50
SX/GLP28	fits 28mm ø screw necks	£5.50
SX/GLP30	fits 30mm ø screw necks	£9.20
SX/GLP34	fits 34mm ø screw necks	£11.40
SX/GLP45	fits 45mm ø screw necks	£10.00

Adaptors, PTFE, screw neck

SX/GLT36	fits 36mm ø screw necks	£46.00
SX/GLT38	fits 38mm ø screw necks	£46.00

Adaptors, polypropylene, standard taper

SX/NSP18.8	fits taper NS19/27	£17.30
SX/NSP24	fits taper NS24/29	£17.30
SX/NSP29.2	fits taper NS29/34	£17.30
SX/NSP45	fits taper NS45/10	£17.30

Adaptors, PTFE, standard taper


SX/NST18.8	fits taper NS19/26	£44.00
SX/NST24	fits taper NS24/29	£44.00
SX/NST29.2	fits taper NS29/34	£44.00

Feed tube, PTFE

Reference	Length x i.d	To fit	Price
SX/09.511.707	300x5mm	520: 2, 5, 10ml	£8.10
SX/09.521.621	300x7.5mm	521: 20ml	£8.10
SX/09.521.622	300x10mm	521: 50ml	£8.10
SX/09.521.623	375x10mm	521: 100ml	£8.10



Adaptors



Amber Glass and Polythene Bottles, square