


Flasks/Funnels



Glass Volumetric Flasks, Class A

Glass volumetric flasks, Class A

Class A

With one mark and interchangeable "Polystop" stopper. BS1792 (ISO 1042). All in borosilicate glass and individually checked for accuracy.

Reference	Capacity	Tolerance ±	Pack quantity	Price
JP/RB258/A/5	5ml	0.025 ml	5	£25.00
JP/RB258/B/5	10ml	0.025 ml	5	£25.00
JP/RB258/D/5	20ml	0.04 ml	5	£29.00
JP/RB258/E/5	25ml	0.04 ml	5	£29.00
JP/RB258/F/5	50ml	0.06 ml	5	£30.00
JP/RB258/G/5	100ml	0.10 ml	5	£41.00
JP/RB258/H/2	200ml	0.15 ml	2	£23.00
JP/RB258/I/2	250ml	0.15 ml	2	£23.00
JP/RB258/J/2	500ml	0.25 ml	2	£31.00
JP/RB258/K/2	1000ml	0.40 ml	2	£40.00
JP/RB258/L/S	2000ml	0.60 ml	1	£36.00
JP/RB258/M/S	5000ml	2 ml	1	£94.00

FUNNELS

Filter funnels

Short stem

These Pyrex conical filter funnels have an accurate 60° bowl angle for precise fitting of filter paper. Produced with beaded rims.



Short Stem Filter Funnels

Reference	Top dia.	Stem	Price
BY/2140/01M	35mm	35mm	£7.00
BY/2140/02M	45mm	45mm	£7.00
BY/2140/03M	55mm	55mm	£7.60
BY/2140/06M	80mm	80mm	£7.70
BY/2140/10M	100mm	100mm	£9.10
BY/2140/11M	125mm	125mm	£12.00
BY/2140/12M	150mm	150mm	£18.40
BY/2140/14M	200mm	175mm	£31.00

Long stem

These Pyrex conical filter funnels are similar to the short stem models but have a stem approximately 150mm long.



Long Stem Filter Funnels

Reference	Top dia.	Stem	Price
BY/2145/03M	55mm	150mm	£16.00
BY/2145/06M	80mm	150mm	£16.90
BY/2145/10M	100mm	150mm	£20.00