

Coliform Swabcheck

The presence of coliforms is indicated by a change from red to yellow. The more rapid the change the higher the level of coliform bacteria. Supplied in packs of 25.

Reference: WN/10498406 £64.00

E.Coli Swabcheck

Used for the detection of E.Coli on surfaces. Supplied in packs of 25

Reference: WN/10498402 £106.00

Total Count Swabcheck

Used for the non-selective development and enumeration of all aerobic bacteria on surfaces in accordance with HACCP. Supplied in packs of 30

Reference: WN/10498315 £99.00



Microbiological Monitors

Single use, pre-sterilised filtering units with welded fixed membranes and culturing devices. Ideal for monitoring contaminants in liquid samples drawn from raw materials to finished products. Supplied in packs of 50.

Reference	Description	Diameter	Price
WN/10497511	0.20um white/black grid	47mm	£78.00
WN/10497500	0.45um white/black grid	47mm	£78.00
WN/10497501	0.45um white/black grid bagged	47mm	£106.00
WN/10497502	0.45um black/white grid	47mm	£78.00
WN/10497503	0.80um black/white grid	47mm	£78.00
WN/10497603	0.20um white/black grid	56mm	£73.00
WN/10497600	0.45um white/black grid	56mm	£90.00
WN/10497601	0.45um black/white grid	56mm	£90.00
WN/10497602	0.80um black/white grid	56mm	£90.00

