

Microbiology/Mixers



Media

Media

Reference	Description	Price
WN/10496714	EC broth 9ml vial/ Durham tube pk20	£18.80
WN/10496709	EC broth with MUG 9ml vial/ Durham tube pk20	£18.80
WN/10496120	Enterococcus broth 2ml ampoule pk50	£56.00
WN/10496722	Lauryl sulphate/tryptose broth 9ml vial pk20	£18.80
WN/10496124	M-FC broth 2ml ampoule pk50	£38.40
WN/10496102	M-TGE media, 2ml ampoule pk50	£34.90
WN/10496119	Pseudomonas broth 2ml ampoule pk50	£38.40
WN/10496113	Total count media with TTC 2ml ampoule pk50	£37.00

MIXERS

Vortex Mixer

Eliminates tedious hand mixing

For rapid and efficient mixing, the Cyclone is far superior to hand shaking. The Cyclone creates a vortex in the tube which gives the highly effective mixing action.

- Variable speed control, with the advantages of controlled mixing over a wide speed range.
- Operates at slow speed for gentle mixing necessary in many applications, where damage to cell structure is to be avoided.
- May be set for continuous running at the chosen speed or on demand by inserting the tube into the cup.
- Skirted rubber cup prevents liquid entering the mixer casing.
- Will accept vessels up to 25mm diameter.
- Excellent stability is achieved through a combination of weight and non-slip feet.



Vortex Mixer

SPECIFICATIONS

Dimensions (wxdxh):	150 x 175 x 150 mm
Watts:	75W
Voltage:	230V, 50Hz

Reference	Price
NO/CM1-1	£163.00