

## IKA Vortex Mixer MS2

Variable speed vortex mixer for fast, thorough mixing of liquids in test tubes up to 30mm diameter. Can be operated continuously, or intermittently when pressure of the tube on the pad will start the mixer. Also includes soft foam platform for processing larger vessels, eg 100ml flasks.

### SPECIFICATIONS

Speed range:	0 to 2500 rpm
Orbit:	4.5mm
Dimensions (wxdxh):	115x225x60mm
Weight:	3.5kg



IKA Vortex Mixer MS2

## Ordering Information

IK/L002050 Vortex mixer MS2 £191.00

## Vortex Mixers, SA7 & SA8

### Vortex, Mixer, fixed speed, SA7

- Fixed speed 2500rpm
- Robust die-cast construction
- Automatic press start

### Vortex, Mixer variable speed, SA8

- Variable speed control between 200 and 2500rpm
- Intermittent or continuous mode
- Stable low profile body
- Robust die-cast construction

Integral retort rod fixing allows vessels to be secured above the vortex action for long-term mixing. The SA8 can mix a variety of other vessels when used in conjunction with the SA8/1 accessory pack detailed overleaf.

### SPECIFICATIONS

	CN/SA7	CN/SA8
Speed	2500rpm	200 to 2500rpm
Orbit diameter	4.2mm	4.2mm
Touch mode	Yes	Yes
Continuous mode	No	Yes
Dimensions (wxdxh)	135 x 215x 78mm	135 x 215 x 78mm
Weight	3.2kg	3.2kg
Electrical supply	90-240V50/60Hz, 20W	90-240V, 50/60Hz, 20W



Vortex Mixers, SA7 & SA8

## Ordering Information

CN/SA7 Vortex Mixer, fixed speed £149.00  
 CN/SA8 Vortex Mixer, variable speed £180.00  
 CN/SR2 Retort rod, 300mm x 12.5mm diameter £13.90