

## PH ELECTRODES

Electrode Type:

P11 P12 P13 P14 P15 P16 P17 P18 P19 P21 P10 01 02 03



### Sentek Electrodes

Competitively priced electrodes offering fast response, low noise and drift, and reliable performance. All electrodes are supplied with one metre cable and BNC plug (except SK/452-75 which is supplied with 2x4mm plugs).

Reference	Type	Temperature range °C	Reference type	Junction type	Dimensions lengthxdia	Body type	Price
<b>Combination electrodes, pH 0-14</b>							
SK/111-75	P11	0-80	AgCl	annular ceramic	120x12mm	glass	£35.00
SK/112-75	P12	0-80	AgCl	frit ceramic	150x6mm	glass	£50.00
SK/113-75	P13	0-80	AgCl	frit ceramic	90x4.5mm	glass	£52.00
SK/114-75	P14	0-60	AgCl	porous teflon	120x12mm	epoxy	£35.00
SK/115-75	P15	0-50	HgCl	annular ceramic	120x12mm	glass	£76.00
SK/116-75	P16	0-50	HgCl	frit ceramic	90x6mm	glass	£82.00
SK/117-75	P17	0-80	AgCl	porous teflon	120x12mm	glass	£82.00
SK/118-75	P18	0-80	AgCl	annular ceramic	120x12mm	glass	£52.00
SK/119-75	P19	0-80	AgCl	frit ceramic	40x6mm	glass	£60.00
SK/121-75	P21	0-80	AgCl	sleeve	120x12mm	glass	£78.00
<b>Glass electrode, pH 0-14</b>							
SK/110-75	P10	0-80	-	-	120x12mm	glass	£26.00
<b>Metal electrodes</b>							
SK/451-75	01	0-80	AgCl	annular ceramic	120x12mm	glass	£66.00
SK/452-75	02	0-80	-	-	120x12mm	glass	£74.00
SK/453-75	03	0-80	-	-	120x12mm	glass	£59.00