

PIPETTES

VOLAC Bulb Pipettes

Class B

Featuring bold figures in permanent amber stain, durable colour codes and "Tufftip" tips to resist damage. All pipettes are calibrated to deliver. To BS1583 (ISO 648).

Reference	Capacity	Accuracy ±	Colour code	Pack quantity	Pack price
JP/R370/A/5	1ml	0.015ml	Blue	5	£10.80
JP/R370/B/5	2ml	0.02ml	Orange	5	£10.80
JP/R370/C/5	3ml	0.03ml	Black	5	£13.80
JP/R370/D/5	4ml	0.03ml	2 red	5	£13.80
JP/R370/E/5	5ml	0.03ml	White	5	£11.80
JP/R370/F/5	10ml	0.04ml	Red	5	£13.80
JP/R370/G/5	15ml	0.05ml	Green	5	£21.00
JP/R370/H/5	20ml	0.06ml	Yellow	5	£20.00
JP/R370/I/5	25ml	0.06ml	Blue	5	£19.70
JP/R370/J/5	50ml	0.1ml	Red	5	£25.00
JP/R370/K/2	100ml	0.15ml	Yellow	2	£14.70

Class A

Guaranteed to give both volumes and flow times within the Class A limits. Suffix NC denotes verification only, suffix WAC denotes pipettes with a Works Accuracy Certificate. All NC and WAC Class A pipettes are individually calibrated and tested, numbered, fully traceable and marked with the time of flow to deliver. To BS1583 (ISO 648).

Reference	Capacity	Accuracy ±	Colour code	Pack quantity	Pack price
JP/R371/NC/A/5	1ml	0.008ml	Blue	5	£39.00
JP/R371/NC/B/5	2ml	0.01ml	Orange	5	£39.00
JP/R371/NC/C/5	3ml	0.015ml	Black	5	£43.00
JP/R371/NC/D/5	4ml	0.015ml	2 Red	5	£43.00
JP/R371/NC/E/5	5ml	0.015ml	White	5	£43.00
JP/R371/NC/F/5	10ml	0.02ml	Red	5	£46.00
JP/R371/NC/G/5	15ml	0.025ml	Green	5	£49.00
JP/R371/NC/H/5	20ml	0.03ml	Yellow	5	£60.00
JP/R371/NC/I/5	25ml	0.03ml	Blue	5	£66.00
JP/R371/NC/J/5	50ml	0.05ml	Red	5	£68.00
JP/R371/NC/K/2	100ml	0.08ml	Yellow	2	£29.00



Volac Bulb Pipettes