


PIPETTE FILLERS

Fast Release Pipette Pumps

Low cost, fast release Pipette Pumps designed for rapid, effective, single-handed operation, with either glass or plastic pipettes. The soft elastic chuck securely holds the pipette in place. Rotation of the thumb wheel draws up liquid, to empty rapidly, apply light thumb pressure to the fast-release lever; for gradual dispensing, rotate the thumb wheel.

Ordering Information

PPD.100/2	Up to 2ml blue	£14.60
PPD.100/10	Up to 10ml green	£14.80
PPD.100/25	Up to 25ml red	£19.80



Fast Release Pipette Pumps

Pipette Controller

Safe, easy pipetting

This controller is the latest in pipetting technology. For fixed volume and graduated pipettes from 0.1 to 100ml capacity, its special valve system provides extremely simple and easy control.

- Suction bellows – simple squeeze action before pipetting, with no need to operate a valve
- Pipetting lever – just one small lever for filling, adjusting the volume, and delivery
- Blowball – to blow the last drop out of blow-out pipettes
- Pipette adapter – removable and internally tapered for secure gripping of all pipettes from 0.1 to 100ml capacity. It has a built-in hydrophobic membrane filter for protection against liquid entering the controller. The adapter and filter are both sterilisable at 121°C

The unit weighs just 95g and is supplied with full instructions.

Ordering Information

Reference	Description	Price
BY/PFL001N	Pipette controller	£46.00
BY/PFL002	Replacement filter	£4.60



Pipette Controller

PVC teats

Approximately 2ml capacity.

Ordering Information

LHW.TS1	Black, pack of 25	£10.80
LHW.TS2	Red, pack of 25	£10.80



PVC teats