

PIPETTORS

Acura® 825/835 Adjustable Volume Micro and Macropipettes

Equipped with a precision digital display, the new line is continuously adjustable down to a very fine volume resolution. Single-handed volume setting is easy and measuring unknown volumes by aspiration is also possible with a high degree of accuracy.

A total of 8 micropipettes covering the range from 0.1 to 1000µL and 2 macropipettes 0.5 - 5ml and 1 - 10ml are available.

Macropipettes - maxi advantages

The 835 models accommodate nozzle filters for added protection against liquid fill and contamination. Optional adapter-nozzle allows fitting of Pasteur pipettes on 5ml model, in addition to Polypropylene tips.

Minimum maintenance/sterilisation

Shock, heat, chemical and UV light resistance is provided by carefully selected materials and good construction. Cleaning and sterilisation can be achieved by autoclaving the instrument, fully assembled, at 121°C for 20 minutes.

Colour coding

The Acura® pipettes have colour-coded display windows and button caps denoting the tip colour. If you require other colour combination for individual, lab or department identification then choose from 14 colours of 'smarties', or pick the 'smartie' mix which contains all colours.


Ordering information

Acura® 825, adjustable volume micropipettes


Reference	Volume µl	Division µl	Tip size	Price
SX/825.0002	0.1-2	0.002	ultra 10µl	£109.00
SX/825.0010	0.5-10	0.01	ultra 10µl	£109.00
SX/825.0010Y	1-10	0.01	200µl	£102.00
SX/825.0020	2-20	0.02	200µl	£102.00
SX/825.0050	5-50	0.1	200µl	£102.00
SX/825.0100	10-100	0.1	200µl	£102.00
SX/825.0200	20-200	0.2	200µl	£102.00
SX/825.1000	100-1000	1	1000µl	£102.00

Disposable tips for SX/815 and SX/825 micropipettes

Reference	Description	Pack size	Price
SX/308A	Standard yellow tips	Pack of 1000	£16.60
SX/309A	Standard blue tips	Pack of 1000	£13.90
SX/302	Micro tips, natural, 10µl	Pack of 100	£3.80
QP/114	Micro tips, natural, 10µl	Pack of 1000	£15.30



Acura® 825/835 Adjustable
Volume Micro and Macropipettes



Acura® 825/835 Adjustable
Volume Micro and Macropipettes



Acura® 825/835 Adjustable
Volume Micro and Macropipettes