

MonoDirect Magnetic Stirrer

A directly controlled electromagnetic stirrer suitable for most general laboratory applications

- Wear-free inductive drive
- Speed range 130 to 1000rpm
- Super flat, only 35mm high
- Stirring volume up to 3 litres
- Safe acceleration to required stirring speed
- Powder-coated stainless steel housing
- Power supply with pre-selectable stirring power
- Illuminated power on lamp

Ordering Information

DH/30110 MonoDirect magnetic stirrer £145.00



MonoDirect Magnetic Stirrer

MaxiDirect Magnetic Stirrer

Electromagnetic stirrer with push button controls and digital display

- Wear-free inductive drive
- Speed range 80 to 2000rpm
- Super flat, only 35mm high
- Stirring volume up to 5 litres
- Safe acceleration to required stirring speed
- Digital display of speed
- 4 power settings, ideal for larger volumes and more viscous liquids
- Memory function stores last setting
- Hermetically sealed stainless steel housing
- Can be cleaned under running water
- Push button controls

Ordering Information

DH/30125 MaxiDirect magnetic stirrer £360.00



MaxiDirect Magnetic Stirrer

SPECIFICATIONS

	MonoDirect	MaxiDirect
Stirring points:	1	1
Stirring volume:	1ml to 3 litres	1ml to 5 litres
Stirring power:	1/5 W	5/10/15/20 W
Speed range:	130 to 1000rpm	80 to 2000rpm
Dimensions:	180x215x35mm	180x215x35mm
Weight:	1.4 kg	2.5 kg
Ambient temperature:	-10 to +40°C at 95% relative humidity	-10 to +40°C at 95% relative humidity
Operating voltage:	8 to 12 V DC	24 V DC
Power supply unit:	power supply unit included	power supply unit included