

# Whatman™ folded filter papers



## Ready-to-use paper solutions

Pre-folded filter papers can save valuable time in the lab, providing ease of use when undertaking repetitive or multiple analyses.

## Customized formats

Whatman™ qualitative and quantitative grade filter papers are now available in new folded cone, pyramid, and flat quadrant formats, in our full range of diameters and grades to support your application needs. There are various stacking and packaging options available to suit customer needs.

## Filter papers: pyramid folded and cone grades

### Specifications

Grade	Nominal thickness (µm)	Nominal basis weight (g/m <sup>2</sup> )	Nominal ash content (%) <sup>1</sup>
1	180	87	0.06
4	210	92	0.06
6	180	100	0.15
40	210	95	0.007
41	215	85	0.007
42	200	100	0.007
540	160	85	0.005

<sup>1</sup> Ash is determined by ignition of the cellulose filter at 900°C in air.



Cone folded filter papers

## Filter papers: cone folded grades

### Ordering information

Diameter (mm)	Grade	Catalog number	Quantity/pack
110	1		1000
110	40		1000
125	40		1000

The above table is an example of products set-up; please contact us at [scientific.support@cytiva.com](mailto:scientific.support@cytiva.com) for a full listing of available products.

Customizable and available in different grades upon request.




Cone folded filter papers

## Filter papers: pyramid folded grades

### Ordering information

Diameter (mm)	Grade	Catalog number	Quantity/pack
90	42		1000
90	540		1000
110	4		1000
110	41		1000
110	42		1000
110	540		1000
125	1		1000
125	4		1000
125	6		1000
125	40		1000
125	41		1000
125	42		1000

Customizable and available in different grades upon request.



Pyramid folded filter papers

## Filter papers: quadrant folded

### Specifications

Grade	Nominal thickness (µm)	Nominal basis weight (g/m <sup>2</sup> )	Nominal ash content (%) <sup>1</sup>
1 FF	180	87	0.06
40 FF	210	95	0.007
41 FF	215	85	0.007
0858 FF	170	75	–

<sup>1</sup> Ash is determined by ignition of the cellulose filter at 900°C in air



Quadrant folded filter papers

## Filter papers: quadrant folded

### Ordering information

Diameter (mm)	Description	Catalog number	Quantity/pack
110	Grade 1 FF Quadrant		500
125	Grade 1 FF Quadrant		500
150	Grade 1 FF Quadrant		500
110	Grade 40 FF Quadrant		500
125	Grade 40 FF Quadrant		500
150	Grade 40 FF Quadrant		500
110	Grade 41 FF Quadrant		500
125	Grade 41 FF Quadrant		500
150	Grade 41 FF Quadrant		500
185	Grade 0858 FF Quadrant		100

Customizable and available in different grades upon request.



44 FF quadrant flat folded filter papers

We also offer pre-pleated filter papers. For more information, please visit:

[cytivalifesciences.com/whatman-RTU-filter-papers](https://cytivalifesciences.com/whatman-RTU-filter-papers)

### cytiva.com

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

Whatman is a trademark of Global Life Sciences Solutions USA LLC or an affiliate.

© 2021 Cytiva

For local office contact information, visit [cytiva.com/contact](https://cytiva.com/contact)

CY23931-13Sep21-FL

