

Thermal Cycler Compatibility

Status as of: 0909

One glance at the table will tell you which BRAND PCR plates are compatible with your thermal cycler. We are continuously updating the table with information from manufacturers and feedback from our customers. Ask for a free sample of our PCR plates to check compatibility with your thermal cycler (www.brand.de), with no obligation. A short feedback from you about the results will help us to complete the table.

	24-, 48-, 96-well, no skirt 7814 11, 7814 15, 7813 50	96-well, elevated skirt 7813 52	96-well, full skirt 7813 53	96-well, half skirt 7814 00	384-well, full skirt 7813 45	384-well, full skirt 7813 47	384-well, full skirt, rigid 7813 48
Applied Biosystems							
2700	●	●		●	●	●	●
3100	●	●				●	●
3130	-	-				●	●
3700	●	●			●	●	●
3730/3730x	●	●				●	●
qPCR 5700		●					
qPCR 7000	●	●					
qPCR 7300	●	●					
qPCR 7500	●	●					
qPCR 7700	●	●					
qPCR 7900 HT		●				●	
9600	●	●		●			
9700	●	●		●	●	●	●
Amersham Biosciences							
MegaBACE 500			●				
MegaBACE 1000			●				
MegaBACE 4000					●	●	●
Biometra							
Uno	●	●	●	●			
Uno II	●	●	●	●	●	●	●
T1 Thermal Cycler	●	●	●	●	●	●	●
T3 Thermal Cycler	-			-			
Tgradient	●	●	●	●			
Trobot	●		●	●	●	●	●
BioRad							
iCycler	●	●	●				
MyCycler	●						
qPCR MyiQ	●		●				
qPCR iQ5	●		●				
Corbett Research							
PalmCycler 96			●				
PalmCycler 384						●	●
Eppendorf							
Mastercycler Gradient	●	●	●	●			
Mastercycler ep	●	●	●				
Mastercycler M 384					●	●	●
Mastercycler ep Realplex			●				

	24-, 48-, 96-well, no skirt 7814 11, 7814 15, 7813 50	96-well, elevated skirt 7813 52	96-well, full skirt 7813 53	96-well, half skirt 7814 00	384-well, full skirt 7813 45	384-well, full skirt 7813 47	384-well, full skirt, rigid 7813 48
Ericomp							
Single Block	●	●					
Twin Block	●	●					
Delta Cycler	●	●					
Hybaid							
Multiblock System MBS	●		●		●	●	●
Omnigene	●	●	●	●	●	●	●
Omn-E	●	●	●	●			
PCR Express	●	●	●	●	●	●	●
PCR Sprint	-		●	-			
pxe	●		●		●	●	●
px2	●		●		●	●	●
Touchdown	●	●	●	●	●	●	●
MJ Research							
BaseStation			●				
qPCR Chromo 4			●				
Dyad/Disciple	●		●			●	●
qPCR Opticon			●				
qPCR Opticon 2			●				
PTC-100	●	●	●	●		●	●
PTC-200	●	●	●	●	●	●	●
PTC-225 Tetrad	●	●	●	●	●	●	●
MWG							
Primus 96	●	●	●	●			
Primus 384					●	●	●
Stratagene							
qPCR Mx4000	●	●					
qPCR Mx3000	●		●				
Robocycler	●	●	●	-	●	●	●
TaKaRa							
TP240			●				
TP3000	●		●				
Techne							
TC-412/Flexigene	●	-	●	●	●	●	●
Genius	●	-	●	●	●	●	●
TC-512/Touchgene Gradient	●	-	●	●	●	●	●
TC-3000X	●*	-	-	-	-	-	-
Transgenomic							
Wave System			●				

* compatible with 7814 11 and 7814 15

● = Real Time PCR compatible ● = compatible - = not compatible □ = no info

