

# FRYKA

FRYKA  
Kältetechnik



## Organizing- and storage systems

### Drawer Shelves RS

for standard Cryoboxes.



Construction with extra spacing for optimal, evenly and fast cooling of the complete freezing material.



Both Systems are available for the different FRYKA- freezer type.  
(see backside).

Special custom made shelves on request.

### Stainless Steel Shelves RN

for standard containers.



**See back for technical data and different versions**

**FRYKA**-Kältetechnik GmbH

Ohmstr. 4 / D-73730 Esslingen / phone ++49-711-310599-0 / fax ++49-711-310599-29

[www.fryka.de](http://www.fryka.de) / [info@fryka.de](mailto:info@fryka.de)

## Drawer Shelves for freezing- and storage boxes

Freezer type	Type of shelf	Size per drawer	Number of drawers
		W x L x H mm	Drawers / boxes
Cold box B 30	B 30-RS	300x280x52	4 / 16
Cold box B 35	B 35-RS	420x280x52	5 / 30
Freezer TS ..-100	TS 100-RS	420x420x52	9 / 81
Freezer TUS ..-100	TS 100-RS	420x420x52	9 / 81

If higher freezing- and storage boxes (than 50mm) are used an appropriate number of drawers has to be removed.

## Stainless Steel Shelves for standard containers

Freezer type	Number per unit + type of shelves	Number of shelves splints
Cold box B 30	1x B 30-RN	6
Freezer TT ..-55	2x TT 55-RN	9

## Standard Container stainless steel CN 18/10

Type	Proportions in mm			Volume in Litres
	Length	Width	Height	
PL 20	300	240	01	--
1/20	300	240	20	--
1/10	300	240	40	1,7
1/7	300	240	65	3,3
1/4	300	240	100	5,3
1/3	300	240	150	8,4
1/2	300	240	200	10,9