

# FRYKA

## Temperature Recorder TSG 30/100

Line writer with text print

Standing case with carrying grip

For roles or folding paper

Feed 2 cm / hr.

Record duration 30 days / folding paper block


Temperature range +0 bis -100°C

Input signal 0 bis -1,0 Volt

Suitable for all FRYKA deep temperature devices with 10mV-output

Dimensions: 34 x 28 x 17 cm ( width x length x height)

Electr. Connection 230 V / 50 Hz



**FRYKA**-Kältetechnik GmbH

Ohmstr. 4 / [73730](http://www.fryka.de) Esslingen / phone 0711-310599-0 / fax 0711-310599-29  
[www.fryka.de](http://www.fryka.de) / [info@fryka.de](mailto:info@fryka.de)